

Cutisorb® Ultra




**Super absorbent
wound dressing**



Cutisorb® Ultra – Superior exudate management that

Cutisorb® Ultra offers an innovative 5-layer approach to exudate management. State-of-the-art materials, including super absorbent particles and a unique distribution layer, help provide superior absorption and keep excess moisture away from skin, even under compression. This, combined with a non-sticking, soft wound contact layer, ensures more gentle and less frequent dressing changes for an improved patient quality of life during treatment.




works under pressure



Superior Fluid Management

- Absorbs more than 4x the exudate of other wound dressings*
- Locks exudate away from the wound and surrounding skin even under pressure
- High fluid capacity results in less frequent dressing changes, lower cost of treatment, and allows longer periods of undisturbed healing



1 Cutisorb® Ultra 2 Standard absorption dressing



Indications

For the management of moderate to highly exuding wounds

- pressure ulcers
- arterial ulcers
- diabetic foot ulcers
- fistulas
- venous leg ulcers
- infected wounds under medical supervision

Benefits to the Caregiver and Patient

Caregiver

- Less frequent dressing changes
- Time and cost savings
- Superior fluid management
- Reliable fluid retention, even under compression



Patient

- Fewer, more gentle dressing changes
- Reduced risk of maceration
- Cushioning comfort
- Improved quality of life

* Data on file

Cutisorb® Ultra

Superior exudate management that works under pressure



Retains fluid under compression



PROTECTS against maceration



Highly absorbent



FEWER dressing changes



Locks in bacteria and MMPs¹



promotes **FASTER** wound healing

Ref. No.	Dressing Size	Dressings/ Box	HCPSCS
72639-00	10 x 10 cm 4" x 4"	20	A6251
72639-01	10 x 20 cm 4" x 8"	20	A6252
72639-02	20 x 20 cm 8" x 8"	20	A6252
72639-03	20 x 30 cm 8" x 12"	20	A6253

Suggested HCPSCS Code


HCPSCS Code listed is only suggested code; other codes may be more applicable depending on product use. Codes are subject to change without notice.

¹ Super absorbent particles

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
To order toll-free: BSN medical (+1) 800 552 1157





Wound Phase

-  Necrotic
-  Infected
-  Sloughy
-  Granulating
-  Epithelializing

Wound Depth

-  Shallow
-  Deep
-  Shallow + deep

Exudate

-  Dry to low
-  Low to moderate
-  Moderate to high
-  High to very high

